


## Anti-Kat13A / Src1 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q15788

Item Number: CM0005898



### Introduction

Rabbit monoclonal antibody targeting human KAT13A/SRC1 (nuclear receptor coactivator 1), suitable for Western blot, immunoprecipitation, and ELISA. Used in nuclear receptor signaling and transcription regulation research. UniProt Q15788.

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | KAT13A / SRC1 Rabbit mAb                                                                                               |
| Remarks/Alias    | SRC1; KAT13A; RIP160; F-SRC-1; bHLHe42; bHLHe74; KAT13A / SRC1                                                         |
| Species          | Human                                                                                                                  |
| Gene ID (Human)  | 8648                                                                                                                   |
| Gene ID          | 8648                                                                                                                   |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 400-750 of human KAT13A / SRC1 (Q15788). |
| Source           | Rabbit                                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                                  |
| Application      | WB, IP, ELISA                                                                                                          |
| Cross Reactivity | Human                                                                                                                  |
| SWISS            | Q15788                                                                                                                 |
| Protein Weight   | 157kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                                |