

Anti-Kat1/Hat1 Rabbit Monoclonal Antibody (Wb, Ip, Ihc-P, If/Icc, Elisa) - O14929

Item Number: CM0011236



Introduction

Rabbit monoclonal antibody against human KAT1/HAT1 (SWISS O14929), validated by knockout and suitable for WB, IP, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] KAT1/HAT1 Rabbit mAb
Remarks/Alias	KAT1; KAT1/HAT1
Species	Human
Gene ID	8520
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-400 of human KAT1/HAT1 (NP_003633.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14929
Protein Weight	50 kDa
Shipping	Ice bag