


Anti-Kat2B/Pcaf Rabbit Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - Q92831

Item Number: CM0003983



Introduction

CamelBio offers Anti-KAT2B/PCAF Rabbit Monoclonal Antibody, validated for WB, IP, IHC-P, ELISA. Detects human, mouse, rat KAT2B/PCAF, a histone acetyltransferase involved in transcriptional activation and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	KAT2B/PCAF Rabbit mAb
Remarks/Alias	CAF; PCAF; P/CAF; KAT2B/PCAF
Species	Human
Gene ID	8850
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-450 of human KAT2B/PCAF (NP_003875.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92831
Protein Weight	93 kDa
Shipping	Ice bag