


Anti-Kat2B/Pcaf Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92831

Item Number: CM0012935



Introduction

Rabbit polyclonal antibody targeting human KAT2B/PCAF, suitable for Western blot, IF/ICC, and ELISA. Detects KAT2B/PCAF in human, mouse, and rat. Supports epigenetics and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	KAT2B/PCAF Rabbit pAb
Remarks/Alias	CAF; PCAF; P/CAF; KAT2B/PCAF
Species	Human
GeneID (Human)	8850
GeneID	8850
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-86 of human KAT2B/PCAF (NP_003875.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92831
Protein Weight	93kDa
Shipping	Ice bag