

Anti-Kat5/Tip60 Polyclonal Antibody For Wb, If/Icc, Elisa - Q92993

Item Number: CM0021805



Introduction

High-quality rabbit polyclonal antibody against human KAT5/Tip60, validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat KAT5 with a predicted molecular weight of 59 kDa. Ideal for chromatin and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	KAT5/Tip60 Rabbit pAb
Remarks/Alias	TIP; ESA1; PLIP; TIP60; cPLA2; HTATIP; ZC2HC5; HTATIP1; NEDFASB; KAT5/Tip60
Species	Human
Gene ID	10524
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human KAT5/Tip60 (NP_006379.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92993
Protein Weight	59kDa
Shipping	Ice bag