

Anti-Kat6B Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wyb5

Item Number: CM0019731



Introduction

Rabbit polyclonal antibody against human KAT6B histone acetyltransferase. Validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Useful for transcription and chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	KAT6B Rabbit pAb
Remarks/Alias	qkf; MORF; MOZ2; GTPTS; MYST4; ZC2HC6B; querkopf; KAT6B
Species	Human
Gene ID	23522
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 400-500 of human KAT6B (NP_001243398.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WYB5
Protein Weight	231kDa
Shipping	Ice bag