


Anti-Kat7 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95251

Item Number: CM0027797



Introduction

High-quality rabbit polyclonal antibody against KAT7 (HBO1/MYST2). Validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat KAT7. Recombinant immunogen, 71kDa, O95251. Ideal for chromatin and leukemia research.

[Learn More](#)

Parameter	Value
Product Name	KAT7 Rabbit pAb
Remarks/Alias	HBO1; HBOA; MYST2; ZC2HC7; KAT7
Species	Human
Gene ID (Human)	11143
Gene ID	11143
Immunogen	Recombinant fusion protein corresponding to amino acids 55-193 of human KAT7 (NP_008998.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95251
Protein Weight	71kDa
Shipping	Ice bag