

Anti-Kat8/Myst1/Mof Monoclonal Antibody For If/Icc, Ihc-P, Elisa - Q9H7Z6

Item Number: CM0010259



Introduction

High-quality rabbit monoclonal antibody targeting histone acetyltransferase KAT8/MYST1/MOF (UniProt Q9H7Z6). Suitable for IF/ICC, IF-P, IHC-P, ELISA applications. Detects human and rat. Useful for epigenetic regulation studies.

[Learn More](#)

Parameter	Value
Product Name	KAT8/MYST1/MOF Rabbit mAb
Remarks/Alias	MOF; hMOF; MYST1; LIGOWS; ZC2HC8; KAT8/MYST1/MOF
Species	Human
Gene ID (Human)	84148
Gene ID	84148
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-458 of human KAT8/MYST1/MOF (Q9H7Z6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9H7Z6
Protein Weight	52kDa
Shipping	Ice bag