

Anti-Kat8/Myst1/Mof Rabbit Polyclonal Antibody For Wb, Elisa - Q9H7Z6

Item Number: CM0022020



Introduction

Rabbit polyclonal antibody against human KAT8/MYST1/MOF for Western blot and ELISA. Cross-reacts with mouse and rat. Targets histone acetyltransferase KAT8, a key epigenetic regulator of H4K16 acetylation and chromatin structure.

[Learn More](#)

Parameter	Value
Product Name	KAT8/MYST1/MOF Rabbit pAb
Remarks/Alias	MOF; hMOF; MYST1; LIGOWS; ZC2HC8; KAT8/MYST1/MOF
Species	Human
GeneID (Human)	84148
GeneID	84148
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human KAT8/MYST1/MOF (NP_892003.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H7Z6
Protein Weight	52kDa
Shipping	Ice bag