

Anti-Kat9/Elp3 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9H9T3

Item Number: CM0003499



Introduction

Rabbit monoclonal antibody against human KAT9/Elp3 (ELP3), validated for WB, IHC-P, ELISA, and human/mouse/rat cross-reactivity. Targets the tRNA acetyltransferase subunit of elongator complex, essential for tRNA modification and neurogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	KAT9/Elp3 Rabbit mAb
Remarks/Alias	KAT9; KAT9/Elp3
Species	Human
GeneID (Human)	55140
GeneID	55140
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-547 of human KAT9/Elp3 (Q9H9T3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H9T3
Protein Weight	62kDa
Shipping	Ice bag