

Anti-Katna1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 075449

Item Number: CM0023353



Introduction

Rabbit polyclonal antibody targeting human KATNA1 (katanin p60 ATPase subunit A1). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Suitable for microtubule severing and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	KATNA1 Rabbit pAb
Remarks/Alias	KATNA1
Species	Human
GeneID (Human)	11104
GeneID	11104
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-167 of human KATNA1 (NP_001191005.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075449
Protein Weight	56kDa
Shipping	Ice bag