

Anti-Katnb1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bva0

Item Number: CM0023752

Introduction

Rabbit polyclonal antibody targeting human KATNB1 (p80 katanin), validated in WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for microtubule severing and cell division studies.

[Learn More](#)

Parameter	Value
Product Name	KATNB1 Rabbit pAb
Remarks/Alias	KAT; LIS6; KATNB1
Species	Human
Gene ID (Human)	10300
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human KATNB1 (NP_005877.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9BVA0
Protein Weight	72kDa
Shipping	Ice bag