

## Anti-Kbtbd3 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Nab2

Item Number: CM0013357

### Introduction

Rabbit polyclonal antibody against human KBTBD3, validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for research on BTB/kelch domain proteins.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | KBTBD3 Rabbit pAb                                                                                                   |
| Remarks/Alias    | BKLHD3; KBTBD3                                                                                                      |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 143879                                                                                                              |
| Gene ID          | 143879                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 520-612 of human KBTBD3 (NP_689646.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q8NAB2                                                                                                              |
| Protein Weight   | 70kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |