

Anti-Kbtbd8 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nfy9

Item Number: CM0022229

Introduction

Affinity-purified rabbit polyclonal antibody against human KBTBD8 (TA-KRP). Suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for E3 ubiquitin ligase complex studies and neural crest specification research.

[Learn More](#)

Parameter	Value
Product Name	KBTBD8 Rabbit pAb
Remarks/Alias	TAKRP; TA-KRP; KBTBD8
Species	Human
GeneID (Human)	84541
GeneID	84541
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-280 of human KBTBD8 (NP_115894.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8Nfy9
Protein Weight	69kDa
Shipping	Ice bag