

Anti-Kcna1 Polyclonal Antibody For Wb, If-P, Elisa - Q09470

Item Number: CM0019243

Introduction

Rabbit polyclonal antibody targeting human KCNA1 (Kv1.1), suitable for Western blot, immunofluorescence (paraffin), and ELISA. Cross-reacts with human, mouse, and rat. Involved in neuronal excitability and kidney magnesium homeostasis.

[Learn More](#)

Parameter	Value
Product Name	KCNA1 Rabbit pAb
Remarks/Alias	EA1; MK1; AEMK; HBK1; HUK1; MBK1; RBK1; KV1.1; KCNA1
Species	Human
Gene ID (Human)	3736
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 426-495 of human KCNA1 (NP_000208.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q09470
Protein Weight	56kDa
Shipping	Ice bag