

Anti-Kcna2 Polyclonal Antibody For Wb, If/Icc, Elisa - P16389

Item Number: CM0026641

Introduction

Rabbit polyclonal antibody against human KCNA2 (Kv1.2, potassium voltage-gated channel subfamily A member 2). Validated for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. Useful for ion transport and neuroscience studies.

[Learn More](#)

Parameter	Value
Product Name	KCNA2 Rabbit pAb
Remarks/Alias	HK4; MK2; HBK5; NGK1; RBK2; DEE32; HUKIV; KV1.2; EIEE32; KCNA2
Species	Human
Gene ID (Human)	3737
Gene ID	3737
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human KCNA2 (NP_004965.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16389
Protein Weight	57kDa
Shipping	Ice bag