

Anti-Kcna4 Polyclonal Antibody For Wb And Elisa - P22459

Item Number: CM0024499

Introduction

Rabbit polyclonal antibody targeting human KCNA4 (Kv1.4) voltage-gated potassium channel. Suitable for Western blot and ELISA; cross-reacts with mouse and rat. Ideal for ion channel and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	KCNA4 Rabbit pAb
Remarks/Alias	HK1; HBK4; PCN2; HPCN2; HUKII; KCNA8; KV1.4; KCNA4L; MCIDDS; KCNA4
Species	Human
Gene ID (Human)	3739
Gene ID	3739
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P22459
Protein Weight	73kDa
Shipping	Ice bag