


## Anti-Kcnab2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13303

Item Number: CM0023229



### Introduction

Rabbit polyclonal antibody against human KCNAB2 (Kv-beta-2). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Detects endogenous human KCNAB2 protein (~41 kDa). Ideal for neuroscience and ion channel research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | KCNAB2 Rabbit pAb                                                                                                |
| Remarks/Alias    | AKR6A5, KCNA2B, HKvbeta2, KV-BETA-2, HKvbeta2.1, HKvbeta2.2, KCNAB2                                              |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 8514                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-40 of human KCNAB2 (NP_003627.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IF/ICC, ELISA                                                                                                |
| Cross Reactivity | Mouse, Rat                                                                                                       |
| SWISS            | Q13303                                                                                                           |
| Protein Weight   | 41kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |