

Anti-Kcnab3 Rabbit Polyclonal Antibody For Wb And Elisa - O43448

Item Number: CM0015801

Introduction

Rabbit polyclonal antibody targeting human voltage-gated potassium channel subunit beta-3 (KCNA3). Validated for Western blot and ELISA, with mouse cross-reactivity. Suitable for neuroscience and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	KCNAB3 Rabbit pAb
Remarks/Alias	AKR6A9; KCNA3B; KCNA3.1B; KV-BETA-3; KCNAB3
Species	Human
GeneID	9196
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 285-404 of human KCNAB3 (NP_004723.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O43448
Protein Weight	44kDa
Shipping	Ice bag