


Anti-Kcnd2 Polyclonal Antibody For Wb, Elisa - Q9Nzv8

Item Number: CM0028268



Introduction

Rabbit polyclonal antibody against KCND2, raised with recombinant protein immunogen (aa 501-630). Validated for WB and ELISA; detects human and mouse KCND2. Useful for studying neuronal excitability and potassium channel function.

[Learn More](#)

Parameter	Value
Product Name	KCND2 Rabbit pAb
Remarks/Alias	RK5; KV4.2; KCND2
Species	Human
Gene ID (Human)	3751
Gene ID	3751
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 501-630 of human KCND2 (NP_036413.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NZV8
Protein Weight	71kDa
Shipping	Ice bag