


Anti-Kcne2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y6J6

Item Number: CM0030350



Introduction

High-quality rabbit polyclonal antibody against KCNE2 (MiRP1), validated for WB and ELISA. Recognizes human and mouse KCNE2, a potassium channel regulatory subunit involved in cardiac IKr current. Ideal for cardiovascular and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	KCNE2 Rabbit pAb
Remarks/Alias	LQT5; LQT6; ATFB4; MIRP1; KCNE2
Species	Human
GeneID (human)	9992
GeneID	9992
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KCNE2 (NP_751951.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y6J6
Protein Weight	14kDa
Shipping	Ice bag