

Anti-Kcnh1 Polyclonal Antibody For Wb, Elisa - 095259

Item Number: CM0028629



Introduction

Rabbit polyclonal antibody against human KCNH1/EAG1 (O95259) for WB and ELISA, cross-reacting with mouse and rat. Immunogen covers amino acids 813-962. Suitable for ion channel research.

[Learn More](#)

Specification	Detail
Product Name	KCNH1 Rabbit pAb
Remarks/Alias	EAG; EAG1; ZLS1; hEAG; TMBTS; h-eag; hEAG1; Kv10.1; KCNH1
Species	Human
GeneID (Human)	3756
GeneID	3756
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 813-962 of human KCNH1 (NP_002229.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95259
Protein Weight	111kDa
Shipping	Ice bag