


Anti-Kcnh2 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q12809

Item Number: CM0027999



Introduction

Rabbit polyclonal antibody against human KCNH2 (HERG, Kv11.1). Validated for WB, IF-P, IHC-P, ELISA with cross-reactivity to mouse. Ideal for cardiac and neurological channel research. SWISS: Q12809.

[Learn More](#)

Parameter	Value
Product Name	KCNH2 Rabbit pAb
Remarks/Alias	ERG1; HERG; LQT2; SQT1; ERG-1; H-ERG; HERG1; Kv11.1; KCNH2
Species	Human
Gene ID (Human)	3757
Gene ID	3757
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 850-950 of human KCNH2 (NP_000229.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12809
Protein Weight	127kDa
Shipping	Ice bag