

Anti-Kcnh6 Polyclonal Antibody For Wb, Elisa - Q9H252

Item Number: CM0021673

Introduction

Rabbit polyclonal antibody against human KCNH6 (Kv11.2), a voltage-gated inwardly rectifying potassium channel. Suitable for Western blot and ELISA applications. Ideal for studying ion transport and cardiac/neuronal signaling.

[Learn More](#)

Parameter	Value
Product Name	KCNH6 Rabbit pAb
Remarks/Alias	ERG2; ERG-2; HERG2; Kv11.2; hERG-2; KCNH6
Species	Human
GeneID	81033
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-220 of human KCNH6 (NP_110406.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H252
Protein Weight	110kDa
Shipping	Ice bag