


Anti-Kcnp2 Polyclonal Antibody For Wb And Elisa - Q9Ns61

Item Number: CM0026233



Introduction

Rabbit polyclonal antibody targeting human KCNIP2 (KChIP2), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for investigating A-type potassium channel modulation in heart and brain.

[Learn More](#)

Parameter	Value
Product Name	KCNIP2 Rabbit pAb
Remarks/Alias	KCHIP2; KCNIP2
Species	Human
GeneID (Human)	30819
GeneID	30819
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human KCNIP2 (NP_775284.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NS61
Protein Weight	31kDa
Shipping	Ice bag