

Anti-Kcnj12 Rabbit Polyclonal Antibody For Wb, Elisa - Q14500

Item Number: CM0016534

Introduction

Rabbit polyclonal antibody to human KCNJ12 (Kir2.2), an ATP-sensitive inward rectifier potassium channel. Validated for WB and ELISA; cross-reacts with mouse and rat. Ideal for studying resting membrane potential and neuronal excitability.

[Learn More](#)

Parameter	Value
Product Name	KCNJ12 Rabbit pAb
Remarks/Alias	IRK2; hIRK; IRK-2; hIRK1; KCNJN1; Kir2.2; kcnj12x; hkir2.2x; KCNJ12
Species	Human
GeneID (Human)	3768
GeneID	3768
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 354-433 of human KCNJ12 (NP_066292.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14500
Protein Weight	49kDa
Shipping	Ice bag