

Anti-Kcnj15 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99712

Item Number: CM0016681

Introduction

High-quality anti-KCNJ15 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA applications. Recognizes human, mouse, and rat KCNJ15 (Kir1.3/Kir4.2), an inward rectifier potassium channel. Ideal for ion channel research.

[Learn More](#)

Parameter	Value
Product Name	KCNJ15 Rabbit pAb
Remarks/Alias	IRKK; KIR1.3; KIR4.2; KCNJ15
Species	Human
GeneID (human)	3772
GeneID	3772
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 221-375 of human KCNJ15 (NP_002234.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99712
Protein Weight	43kDa
Shipping	Ice bag