

Anti-Kcnj2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P63252

Item Number: CM0014251



Introduction

Anti-KCNJ2 rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat Kir2.1 inward rectifier potassium channel. Suitable for cardiovascular and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	KCNJ2 Rabbit Polyclonal Antibody
Remarks/Alias	IRK1, LQT7, SQT3, ATFB9, HHIRK1, KIR2.1, HHBIRK1, KCNJ2
Species	Human
GeneID (human)	3759
GeneID	3759
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 318-427 of human KCNJ2 (NP_000882.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63252
Protein Weight	48kDa
Shipping	Ice bag