


Anti-Kcnj4 Polyclonal Antibody For Wb, Ihc-P, Elisa - P48050

Item Number: CM0016396



Introduction

High-quality Anti-KCNJ4 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Targets human, mouse, rat KCNJ4 inward rectifier potassium channel (Kir2.3). Ideal for ion channel research.

[Learn More](#)

Parameter	Value
Product Name	Anti-KCNJ4 Polyclonal Antibody
Remarks/Alias	HIR; HRK1; IRK3; HIRK2; IRK-3; Kir2.3; KCNJ4
Species	Human
GeneID (Human)	3761
GeneID	3761
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human KCNJ4 (NP_004972.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48050
Protein Weight	50kDa
Shipping	Ice bag