

Anti-Kcnj5 Polyclonal Antibody For Wb, Elisa - P48544

Item Number: CM0029244

Introduction

Rabbit polyclonal antibody against human KCNJ5 (GIRK4), suitable for Western blot and ELISA. Reacts with human, mouse, and rat. Target is a G protein-activated inward rectifier potassium channel involved in cardiac and endocrine function. Ideal for cardiovascular and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	KCNJ5 Rabbit pAb
Remarks/Alias	CIR; GIRK4; KATP1; LQT13; KIR3.4; KCNJ5
Species	Human
Gene ID (Human)	3762
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human KCNJ5 (NP_000881.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48544
Protein Weight	48kDa
Shipping	Ice bag