


Anti-Kcnk6 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y257

Item Number: CM0019296



Introduction

High-quality rabbit anti-KCNK6 polyclonal antibody validated for Western blot and ELISA. Reacts with human, mouse, and rat. Suitable for potassium channel function and NLRP3 inflammasome research.

[Learn More](#)

Parameter	Value
Product Name	KCNK6 Rabbit pAb
Remarks/Alias	TOSS; KCNK8; TWIK2; K2p6.1; TWIK-2; KCNK6
Species	Human
Gene ID (Human)	9424
Gene ID	9424
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 254-313 of human KCNK6 (NP_004814.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y257
Protein Weight	34kDa
Shipping	Ice bag