


Anti-Kcnk9 Polyclonal Antibody For Wb And Elisa - Q9Npc2

Item Number: CM0016569



Introduction

High-quality rabbit polyclonal antibody against KCNK9 (TASK-3) for Western blotting and ELISA; detects human and mouse samples, with a predicted weight of 42 kDa; ideal for ion channel and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	KCNK9 Rabbit pAb
Remarks/Alias	KT3.2; TASK3; BIBARS; K2p9.1; TASK-3; TASK32; KCNK9
Species	Human
Gene ID (Human)	51305
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 265-374 of human KCNK9 (NP_001269463.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NPC2
Protein Weight	42kDa
Shipping	Ice bag