

Anti-Kcnn3/Sk3 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Ugi6

Item Number: CM0007501

Introduction

High-specificity rabbit monoclonal antibody targeting human KCNN3/SK3 (KCa2.3), suitable for WB and ELISA; cross-reacts with rat. Predicted MW 81 kDa, UniProt Q9UGI6.

[Learn More](#)

Product Name	KCNN3/SK3 Rabbit mAb
Remarks/Alias	SK3; ZLS3; hSK3; SKCA3; KCa2.3; KCNN3/SK3
Species	Human
Gene ID	3782
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 631-731 of human KCNN3/SK3 (NP_002240.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9UGI6
Protein Weight	81 kDa
Shipping	Ice bag