

Anti-Kcnn4 Rabbit Polyclonal Antibody For Wb, Elisa - O15554

Item Number: CM0024693



Introduction

Rabbit polyclonal antibody against human KCNN4 (KCa3.1/IKCa1), validated for WB and ELISA. Suitable for studying calcium-activated potassium channels in immunology, vascular biology, and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	KCNN4 Rabbit pAb
Remarks/Alias	IK; IK1; SK4; DHS2; KCA4; hSK4; IKCA1; hKCa4; KCa3.1; hIKCa1; KCNN4
Species	Human
GeneID (Human)	3783
GeneID	3783
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15554
Protein Weight	48kDa
Shipping	Ice bag