


Anti-Kcnq1 Rabbit Polyclonal Antibody For Wb And Elisa - P51787

Item Number: CM0023338



Introduction

KCNQ1 Rabbit Polyclonal Antibody validated for WB and ELISA in Human, Mouse, Rat. Targets potassium channel subunit Kv7.1 involved in cardiomyocyte excitability and ion transport. Ideal for cardiac and neurological research.

[Learn More](#)

Parameter	Value
Product Name	KCNQ1 Rabbit pAb
Remarks/Alias	LQT; RWS; WRS; LQT1; SQT2; ATFB1; ATFB3; JLNS1; KCNA8; KCNA9; Kv1.9; Kv7.1; KVLQT1; KCNQ1
Species	Human
Gene ID (Human)	3784
Gene ID	3784
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-549 of human KCNQ1 (NP_861463.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51787
Protein Weight	75kDa
Shipping	Ice bag