

Anti-Kcnq2 Rabbit Polyclonal Antibody For Wb, Elisa - O43526

Item Number: CM0021846



Introduction

KCNQ2 Rabbit Polyclonal Antibody validated for Western Blot (WB) and ELISA. Cross-reactive with human, mouse, and rat. For studying neuronal excitability and the M-current. O43526.

[Learn More](#)

Parameter	Value
Product Name	KCNQ2 Rabbit pAb
Remarks/Alias	EBN, BFNC, DEE7, EBN1, ENB1, HN5PC, KV7.2, KCNA11, KCNQ2
Species	Human
GeneID (Human)	3785
GeneID	3785
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 466-665 of human KCNQ2 (NP_742105.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43526
Protein Weight	96kDa
Shipping	Ice bag