

Anti-Kcns2 (Kv9.2) Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9ULS6

Item Number: CM0014332

Introduction

Rabbit polyclonal antibody targeting KCNS2 (KV9.2), a delayed-rectifier potassium channel regulatory subunit. Validated for WB, IF-P, ELISA; cross-reacts with human, mouse, and rat. Suitable for ion channel and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	KCNS2 Rabbit pAb
Remarks/Alias	KV9.2; KCNS2
Species	Human
Gene ID (Human)	3788
Gene ID	3788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-187 of human KCNS2 (NP_065748.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9ULS6
Protein Weight	54kDa
Shipping	Ice bag