

Anti-Kcns3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bq31

Item Number: CM0027487

Introduction

Research-grade rabbit polyclonal antibody targeting KCNS3 (Kv9.3), validated in WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. For ion channel regulation studies.

[Learn More](#)

Parameter	Value
Product Name	KCNS3 Rabbit pAb
Remarks/Alias	KV9.3; KCNS3
Species	Human
GeneID (Human)	3790
GeneID	3790
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human KCNS3 (NP_002243.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQ31
Protein Weight	56kDa
Shipping	Ice bag