


Anti-Kctd15 Polyclonal Antibody For Wb And Elisa - Q96Si1

Item Number: CM0026591



Introduction

Rabbit polyclonal antibody against human KCTD15. Recombinant fusion protein immunogen; validated for WB and ELISA. Cross-reacts with mouse and rat. Detects ~32 kDa band. Useful for studying neural crest formation and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	KCTD15 Rabbit pAb
Remarks/Alias	KCTD15
Species	Human
Gene ID (Human)	79047
Gene ID	79047
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-234 of human KCTD15 (NP_076981.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96Si1
Protein Weight	32kDa
Shipping	Ice bag