

Anti-Kctd9 Monoclonal Antibody For Wb, Elisa - Q7L273

Item Number: CM0003595



Introduction

Rabbit monoclonal antibody targeting human KCTD9 (BTBD27), a BTB/POZ domain-containing E3 ubiquitin ligase adapter. Validated for WB and ELISA in human samples. Ideal for studying protein degradation pathways.

[Learn More](#)

Parameter	Value
Product Name	KCTD9 Rabbit mAb
Remarks/Alias	BTBD27; KCTD9
Species	Human
Gene ID (Human)	54793
Gene ID	54793
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KCTD9 (Q7L273)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q7L273
Protein Weight	43 kDa
Shipping	Ice bag