


Anti-Kdm1 (Lsd1) Rabbit Polyclonal Antibody For Wb, Ihc-P - O60341

Item Number: CM0024250



Introduction

Rabbit polyclonal antibody against human KDM1 (LSD1) validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for chromatin and epigenetic research.

[Learn More](#)

Parameter	Specification
Product Name	KDM1 Rabbit pAb
Remarks/Alias	AOF2; CPRF; KDM1; LSD1; BHC110
Species	Human
Gene ID (Human)	23028
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-380 of human KDM1 (NP_055828.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60341
Protein Weight	93kDa
Shipping	Ice bag