


Anti-Kdm1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - O60341

Item Number: CM0015589



Introduction

CamelBio offers high-quality Anti-KDM1 Rabbit Polyclonal Antibody (O60341) for WB, IP, ELISA, with broad human, mouse, rat reactivity. Ideal for epigenetic regulation, cancer, and developmental biology studies.

[Learn More](#)

Parameter	Value
Product Name	KDM1 Rabbit pAb
Remarks/Alias	AOF2; CPRF; KDM1; LSD1; BHC110
Species	Human
Gene ID (Human)	23028
Gene ID	23028
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 365-515 of human KDM1 (NP_055828.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60341
Protein Weight	93kDa
Shipping	Ice bag