


Anti-Kdm1A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O60341

Item Number: CM0022749



Introduction

KDM1A rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Specific to human, mouse, rat lysine-specific histone demethylase 1A. Essential for epigenetic regulation and gene expression research. - O60341

[Learn More](#)

Parameter	Value
Product Name	KDM1A Rabbit pAb
Remarks/Alias	AOF2; CPRF; KDM1; LSD1; BHC110; KDM1A
Species	Human
GeneID (Human)	23028
GeneID	23028
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 130-380 of human KDM1A (NP_055828.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60341
Protein Weight	93kDa
Shipping	Ice bag