


Anti-Kdm1/Lsd1 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - O60341

Item Number: CM0006751



Introduction

High-quality rabbit monoclonal antibody against human KDM1/LSD1, validated for Western blot, immunofluorescence, IHC, and ELISA. Reacts with human, mouse, rat, and monkey. Ideal for epigenetics and cancer research.

[Learn More](#)

Parameter	Value
Product Name	KDM1/LSD1 Rabbit mAb
Remarks/Alias	AOF2; CPRF; KDM1; LSD1; BHC110; KDM1/LSD1
Species	Human
GeneID (Human)	23028
GeneID	23028
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-275 of human KDM1/LSD1(NP_055828.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	O60341
Protein Weight	93kDa
Shipping	Ice bag