

Anti-Kdm2A Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y2K7

Item Number: CM0023008



Introduction

Rabbit polyclonal anti-KDM2A antibody, raised against a synthetic peptide from the N-terminus of human KDM2A. Validated for WB and ELISA, cross-reacts with human and mouse. Targets a histone demethylase that regulates chromatin and circadian rhythm.

[Learn More](#)

Parameter	Specification
Product Name	KDM2A Rabbit pAb
Remarks/Alias	FBL7; CXXC8; FBL11; FBXL11; JHDM1A; LILINA; KDM2A
Species	Human
GeneID (human)	22992
GeneID	22992
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KDM2A (NP_036440.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y2K7
Protein Weight	133kDa
Shipping	Ice bag