

Anti-Kdm2B Polyclonal Antibody For Wb And Elisa - Q8Nhm5

Item Number: CM0019728



Introduction

KDM2B Rabbit pAb is a polyclonal antibody for western blot and ELISA, targeting human lysine-specific demethylase 2B. Ideal for epigenetic research and chromatin biology studies.

[Learn More](#)

Parameter	Value
Product Name	KDM2B Rabbit pAb
Remarks/Alias	CXXC2; Fbl10; PCCX2; FBXL10; JHDM1B; KDM2B
Species	Human
GeneID (Human)	84678
GeneID	84678
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 600-700 of human KDM2B (XP_011537177.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8NHM5
Protein Weight	153kDa
Shipping	Ice bag