

Anti-Kdm2B Polyclonal Antibody (Kd Validated) For Wb And Elisa - Q8Nhm5

Item Number: CM0025014



Introduction

Rabbit polyclonal antibody against human KDM2B (Lysine-specific demethylase 2B), validated for Western blot and ELISA. Suitable for chromatin biology and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] KDM2B Rabbit pAb
Remarks/Alias	CXXC2; Fbl10; PCCX2; FBXL10; JHDM1B
Species	Human
Gene ID	84678 (human)
Immunogen	Recombinant fusion protein (aa 1059-1336 of human KDM2B)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8NHM5
Protein Weight	153 kDa
Shipping	Ice bag