


Anti-Kdm3B Polyclonal Antibody For Wb, Ihc-P, Elisa - Q7Lbc6

Item Number: CM0019208



Introduction

Rabbit polyclonal antibody against human KDM3B (Lysine-specific demethylase 3B). Validated for WB, IHC-P, ELISA. Useful for chromatin regulation and epigenetic research. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	KDM3B Rabbit pAb
Remarks/Alias	5qNCA; DIJOS; NET22; CSorf7; JMJD1B; KDM3B
Species	Human
GeneID (Human)	51780
GeneID	51780
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 400-600 of human KDM3B (NP_057688.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q7LBC6
Protein Weight	192kDa
Shipping	Ice bag