

Anti-Kdm4A Monoclonal Antibody For Wb, Elisa - 075164

Item Number: CM0008601



Introduction

Rabbit monoclonal antibody against human KDM4A (JMJD2A), a histone demethylase important in chromatin regulation. Validated for Western blot (WB) and ELISA; cross-reacts with human and mouse. Suitable for epigenetic and gene expression studies.

[Learn More](#)

Parameter	Value
Product Name	KDM4A Rabbit mAb
Remarks/Alias	JMJD2; JHDM3A; JMJD2A; TDRD14A; KDM4A
Species	Human
Gene ID (Human)	9682
Gene ID	9682
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-650 of human KDM4A (075164)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	075164
Protein Weight	121kDa
Shipping	Ice bag