

Anti-Kdm4A Rabbit Monoclonal Antibody For Wb, Ip, Elisa - O75164

Item Number: CM0006727



Introduction

Rabbit monoclonal KDM4A antibody validated for WB, IP, ELISA in human, mouse, rat. Targets histone demethylase JMJD2A, a key epigenetic regulator. Suitable for IVD assay development, chromatin biology studies. 121kDa protein weight.

[Learn More](#)

Parameter	Value
Product Name	KDM4A Rabbit mAb
Remarks/Alias	JMJD2; JHDM3A; JMJD2A; TDRD14A; KDM4A
Species	Human
Gene ID (Human)	9682
Gene ID	9682
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 382-459 of human KDM4A (NP_055478.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75164
Protein Weight	121kDa
Shipping	Ice bag